

HEMORRHOIDAL BAND LIGATOR

Used to ligate internal hemorrhoids, facilitated by an anoscope. Supplied unsterile. For single use only.



Key features:

- Anoscope with in-built LED light
- Anoscope available in two sizes – Medium and Large
- Supplied with 4 latex free rubber bands
- Suction control at finger tip
- Automatic Band counter

Classification:

CATALOGUE
NUMBERS

DEVICE DESCRIPTION

AHBL-04

ACE Hemorrhoidal Band Ligator 04 Bands Without Anoscope

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com acemedicaldevices.com

AHBL-04AM ACE Hemorrhoidal Band Ligator 04 Bands With Medium Anoscope
AHBL-04AL ACE Hemorrhoidal Band Ligator 04 Bands With Large Anoscope