


HEMOCLIP

Hemoclip is intended for endoscopic clip placement within the gastrointestinal tract for the purpose of haemostasis for mucosal / submucosal defects, bleeding ulcers, arteries and polyps.



Key features:

- Efficient, smooth and easy to use handle
- 360° Rotatable clip design
- Atraumatic tip prevents damage to the endoscope
- Unmatched repositionability and improved visualisation
- Opening width 135

Specifications

Size

Diameter 15 mm

Length 2300 mm

Tip Length 8~16 mm