


MECHANICAL LITHOTRIPTOR

Used for regular or emergency Lithotripsy to crush C.B.D. stones endoscopically. Reusable device supplied non sterile.



Key features:

- Available as TTS and OTW variants
- Spirals ensure smooth and easy lithotripsy
- Basket shape hexagonal
- Spiral OTW : 4 mm, length 100 cms
- Spiral TTS : 3.5 mm, length 195 cms
- Reusable design that can be steam sterilized
- Supplied in a carrying case

Specifications

Size

Basket Diameter 25, 30 and 35 mm

Length 400 cm

