

# **NASAL BILIARY DRAINAGE CATHETER (NBDC)**

Used for temporary drainage of the Biliary duct through the use of an indwelling catheter. Supplied sterile in a peel open pouch. Designed for single use only.



## **Key features:**

- Pigtail shape for better retention.
- Each set includes nasal transfer tube and a connector to facilitate irrigation and drainage
- Available as Straight, Straight with Pigtail and Angled with Pigtail variants

## **Classification:**

## Specifications

### Size

**Diameter** 5, 6, 7, 8, 8.5 & 10

Fr

**Length** 250 cm