


NEPHROSTOMY BALLOON DILATOR

Used to dilate the musculofascial tract and renal capsule establishing a Nephrostomy tract during a percutaneous procedure. Designed for single use only.



Key features:

- Exceptional dilation capacity
- No kinking and bending during insertion, and no risk of trauma.
- Fast deflation
- Short, flexible, radio-opaque tip

Classification:

30 Fr Amplatz Sheath provided

Shaft Size: 7 Fr

Balloon Width: 10 mm

Balloon length: 10 cm

Shaft Length : 60 cm