


## PCN Catheter- Pigtail with Suture

Used for percutaneous placement of a Pigtail Catheter in the renal pelvis for Nephrostomy drainage. The locking mechanism prevents the accidental slippage of the catheter. Designed for single use only.



### Key features:

- Highly radiopaque for better visualisation
- The pigtail tip configuration for the retention of the distal tip in the renal pelvis
- Locking type has more retention and can be kept for a longer duration
- Supplied with / without hydrophilic coating

### Specifications

#### Size

**Diameter** 6, 7, 8, 8.5, 9, 10, 12, 14, 16 and 18  
Fr

**Length** 20 - 40 cm