


## PCN CATHETER - MALECOT

Used for percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage. Designed for single use only.



### Key features:

- Highly radiopaque for better visualization
- The Malecot flower design enhances drainage and promotes catheter retention
- Round tip malecot catheters reduces the trauma to the tissues
- Highly suited for drainage of thick viscous fluids
- Supplied with / without hydrophilic coating

### Classification:

Normal Malecot

Suprapubic Malecot

Malecot Catheter with detachable hubs.

## Specifications

### Size

**Diameter** 6, 8, 10, 12, 14, 16, 18, 20 and 22 Fr

**Length** 20 - 40 cm