




## PCN CATHETER – PIGTAIL

Used for the temporary or permanent drainage and irrigation of urine from the kidney in presence of distal obstruction. Designed for single use only.



### Key features:

- Highly radiopaque for better visualization.
- Pigtail tip configuration for the retention of the distal tip in the renal pelvis.
- Soft and flexible material for patient's comfort.
- Supplied with / without hydrophilic coating.

### Classification:

Available with / without Hydrophilic coating

### Specifications

## Size

**Diameter** 6, 7, 8, 8.5, 9, 10, 12, 14, 16, 18 Fr

**Length** 20 - 40 cm