


## PCN DIRECT PUNCTURE SET - MALECOT

Used for percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage. Designed for single use only.



### Key features:

- Highly radiopaque for better visualization.
- The Malecot flower design enhanced drainage and promotes catheter retention.
- Round tip malecot catheter reduces the trauma to the tissue.

### Classification:

Malecot Catheter is supplied along with a Puncture Needle.

Available with / without Hydrophilic Coating.

## Specifications

### Size

**Diameter** 8,9, 10, 12, 14, 16, 18, 20 Fr

**Length** 22 and 30 cm