


## PCN SET - MALECOT

Used for Percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage. Designed for single use only.



### Key features:

- Highly radiopaque for better visualization.
- The Malecot flower design enhances drainage and promotes catheter retention.
- Round tip malecot catheters reduces the trauma to the tissues.
- Highly suited for drainage of thick viscous fluids.
- Malecot Catheter can be supplied with / without Hydrophilic coating

### Classification:

PCN Malecot Set consists of:

- Malecot Catheter
- Set of Fascial Dilators
- Initial Puncture Needle
- "J" Tip Guidewire

- Urine Bag Connector
- Scalpel

Available in Tray packaging as well as pouch packaging.

## Specifications

### Size

**Diameter** 6, 8, 10, 12, 14, 16, 18, 20,22 Fr

**Length** 20 - 35 cm