


PCN SET - PIGTAIL

PCN set is a complete set for percutaneous placement of a Pigtail Drainage Catheter in the renal pelvis for Nephrostomy drainage. Designed for single use only



Key features:

- Highly radiopaque for better visualization
- Pigtail tip configuration for the retention of the distal tip in the renal pelvis
- Looking type has more retention and can kept for a longer duration
- Soft and flexible material for patient's comfort
- Available as PCN catheter with Lock
- Pigtail Catheter can be supplied with / without Hydrophilic coating

Classification:

PCN Pigtail Set consists of:

- Pigtail Catheter
- Set of Fascial Dilators
- Initial Puncture Needle
- "J" Tip Guidewire
- Urine Bag Connector
- Scalpel

Available in Tray packaging as well as Pouch Packaging

Specifications

Size

Diameter 6, 7, 8, 8.5, 9, 10, 12, 14, 16, 18 Fr

Length 22 and 30 cm