


## PEG KIT

Used for Percutaneous Endoscopic Gastrostomy (PEG) placement to provide enteral nutrition to the patient requiring nutritional support. Set consists of Pull PEG tube, bolster, Cold Snare, Looped Guidewire, Silicone connectors, introducer needle with sheath, drape, gauze, clamp and scalpel. This device is supplied sterile in a blister pack. Intended for single use only.



### Key features:

- Kink resistant medical grade Silicon tube
- Safe, innovative nutritional delivery device
- Specially designed Bolster to prevent migration
- ENFit compliant device - ISO-80369-3 is the international standard created to establish secure and unique designs for enteral medical devices and accessories

### Specifications

#### Size

**Diameter** 16, 20, 24 and 28 Fr

**Length** 600 mm