




PLASMA BUTTON - HEMI ELECTRODE

Ace Hemi Electrode is the Plasma Button Electrode used for bipolar vaporization for treatment of benign prostate hyperplasia (BPH). Product supplied Sterile and single use only.



Key features:

- Excellent hemostasis.
- High quality insulation.
- Proprietary design (Patent pending).
- Improved visibility and enhanced durability during procedures.
- Available in two shapes - Round and Oval shape

Classification:

Compatible with:

STORZ 1 Stem Instrument

Olympus 1 Stem Instrument

Wolf 1 Stem Instrument

Other properties:

Available with / without bipolar cable.

Specifications

Size

Diameter 24 Fr