

# PLASMA BUTTON - HEMI ELECTRODE

Ace Hemi Electrode is the Plasma Button Electrode used for bipolar vaporization for treatment of benign prostate hyperplasia (BPH). Product supplied Sterile and single use only.



## Key features:

- Excellent hemostasis.
- High quality insulation.
- Proprietary design (Patent pending).
- Improved visibility and enhanced durability during procedures.
- Available in two shapes - Round and Oval shape

## Classification:

Compatible with:

STORZ 1 Stem Instrument

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email [info@acemedicaldevices.com](mailto:info@acemedicaldevices.com) [acemedicaldevices.com](http://acemedicaldevices.com)

Olympus 1 Stem Instrument

Wolf 1 Stem Instrument

## **Other properties:**

Available with / without bipolar cable.

## **Specifications**

*Size*

**Diameter** 24

Fr