


PLASTIC BILIARY STENTS

Ace Biliary Stents are used to drain obstructed Biliary ducts. This device is supplied sterile in a peel-open pouch. Intended for single use only.



Key features:

- Stents made of polyethylene and / or PTFE
- Smooth cannulation for less trauma to the tissue.
- Highly Radiopaque for better visualization.
- Rounded Taper ends.
- Flap and Pigtail design for better retention.
- High Flow Rate

Classification:

- Standard
- Amsterdam
- Straight

- Tennabaum
- Double Pigtail
- Single Pigtail
- Pancreatica

Specifications

Size

Diameter 5, 7, 8.5, 9,10, 11.5 & 12 Fr

Length 3-20 cm