


PTFE GASTRO GUIDEWIRE



Key features:

- PTFE coated guidewire with 5.5 cm hydrophilic coated - for atraumatic advancement through gastrointestinal tract
- Nitinol with PTFE jacket - easy tracking of stents and instruments
- Excellent torqueability for precise tip control
- Better pushability for increased placement control
- Extremely smooth.
- Safe vessel access.

Specifications

Size

Diameter 0.018, 0.021, 0.025, 0.028, 0.032, 0.035 Inch

Length 260 and 450 cm