


RE - ENTRY MALECOT CATHETER

Used as a combination of nephrostomy catheter with stent for drainage, tract maturation and maintenance for percutaneous procedures.



Key features:

- Highly radiopaque for better visualization
- smooth surfaces for ease of introduction
- uniform taper to reduce trauma.

Specifications

Size

Diameter 8/4, 10/4, 10/5, 10/6, 16/6, 16/7, 12/6, 14/6, 14/7, 16/7, 16/8, 18/8, 20/8 Fr

Length 30 - 48 cm