




ROLLER ELECTRODE - BIPOLAR

Bipolar Roller Electrode is used to vaporize the prostate tissue using high power during TURP Surgery for treatment of BPH. Product supplied Sterile and single use only.



Key features:

- High quality insulation.
- Specially designed for use with bipolar current.
- Compatible with standard resectoscopes.

Classification:

Compatible with

- STORZ 2 Stem Instrument
- Wolf 1 Stem Instrument

Other properties:

Available with / without bipolar cable.

Specifications

Size

Diameter 24 Fr