


STONE EXTRACTION BALLOON

Used for endoscopic extraction of stones in the biliary system and for contrast injection.
Designed for single use only.



Key features:

- High quality balloons designed for safe and efficient procedure
- Triple balloon inflation diameter in a single procedure for dilated ducts
- Injection point above or below the balloon
- Radiopaque bands to identify balloon location fluoroscopically
- Reference marks on the syringe for each of the balloon inflation diameters
- Atraumatic tip prevents damage to the endoscope
- Accepts 0.035" Ace guidewire

Specifications

Size

Diameter 9-12-15 mm

Sheath Diameter 7 Fr

Balloon Inflated Diameter 16 mm

Working Length 2000 mm