

# STONE EXTRACTION BALLOON

Used for endoscopic extraction of stones in the biliary system and for contrast injection. Designed for single use only.



## Key features:

- High quality balloons designed for safe and efficient procedure
- Triple balloon inflation diameter in a single procedure for dilated ducts
- Injection point above or below the balloon
- Radiopaque bands to identify balloon location fluoroscopically
- Reference marks on the syringe for each of the balloon inflation diameters
- Atraumatic tip prevents damage to the endoscope
- Accepts 0.035" Ace guidewire

## Specifications

Size

**Diameter** 9-12-15 mm

**Sheath Diameter** 7 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email [info@acomedicaldevices.com](mailto:info@acomedicaldevices.com) [acomedicaldevices.com](http://acomedicaldevices.com)

**Balloon Inflated Diameter** 16  
mm

**Working Length** 2000 mm