

STRIPED GASTRO GUIDEWIRE

Used to facilitate the placement of instruments during diagnostic or treatment procedures. Designed for single use only.



Key features:

- Straight, Flexible tip - for atraumatic advancement through gastrointestinal strictures
- Nitinol core with striped jacket - easy tracking of stents and instruments
- Better pushability for increased placement control
- Extremely smooth.
- Safe vessel access.
- Available in Straight and Angled Tip

Specifications

Size

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com acomedicaldevices.com

Diameter 0.018, 0.021, 0.025, 0.028, 0.032, 0.035
Inch

Length 260 and 450 cm