

SUPERNOMA BAND ELECTRODE

Bipolar Band Electrode is used to cut and coagulate prostate tissue during TURP surgery. Product supplied Sterile and single use only.



Key features:

- Unique and effective design
- High quality insulation.
- Specially designed for use with bipolar current.
- Compatible with all standard resectoscopes.

Classification:

Compatible with:

- Olympus 1 Stem Instrument

Other properties:

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@ecomedicaldevices.com ecomedicaldevices.com

Available with / without bipolar cable.

Specifications

Size

Diameter 24
Fr