


## SUPRAPUBIC MALECOT CATHETER (SPC)

Used to provide nephrostomy drainage by percutaneous placement of a malecot catheter in the renal pelvis.



### Key features:

- Highly radiopaque for better visualization
- The Malecot flower design enhanced drainage & promote catheter retention
- Round tip malecot catheters reduces the trauma to the tissue

### Classification:

Malecot Catheter supplied with a metal trocar.

### Specifications

#### Size

**Diameter** 6, 8, 10, 12, 14, 16, 18, 20, 22 Fr

**Length** 20 - 40 cm