


CYSTO CATH - SUPRAPUBIC DRAINAGE SET

Ace Cysto Cath is used for percutaneous bladder drainage using a balloon catheter. Intended for single use only.



Key features:

Allows easy and comfortable catheter insertion.

Minimum tissue damage and trauma.

Classification:

- Suprapubic Trocar with peel away sheath or with metallic split cannula.
- Scalpel
- Foley Balloon Catheter (available in Silicone and latex)

Other properties:

Packaging:

In Tray packaging / Pouch packaging

Specifications

Size

Diameter 8, 10, 12, 14, 16, 18 Fr

Length 25 and 40 cm