




## **CYSTO CATH - SUPRAPUBIC DRAINAGE SET**

Ace Cysto Cath is used for percutaneous bladder drainage using a balloon catheter. Intended for single use only.



### **Key features:**

Allows easy and comfortable catheter insertion.

Minimum tissue damage and trauma.

### **Classification:**

- Suprapubic Trocar with peel away sheath or with metallic split cannula.
- Scalpel
- Foley Balloon Catheter (available in Silicone and latex )

## Other properties:

### Packaging:

In Tray packaging / Pouch packaging

## Specifications

### Size

**Diameter** 8, 10, 12, 14, 16, 18 Fr

**Length** 25 and 40 cm