

STONE BASKET - TIPLESS

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for single use only.



Key features:

- Made of kink resistant Nitinol wire
- Handle design allows for complete control.
- Handle disassembly possible
- Basket sheath can be provided in PTFE, PEEK, Polyamide and Braided shaft to suit budget

Specifications

Size

Diameter 1.9, 2.2, 2.4 and 3

Length 70,90,110,120 and 130
cm

