


## STONE BASKETS - SUPER CLENCH NITINOL RETRIEVER

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



### Key features:

- Multiwire design for tight clenching of fragments
- Partial closure forms a web of wires preventing retropulsion of stones
- Sturdy and kink resistant
- Unique design

### Specifications

#### Size

**Diameter** 1.9,2.2,2.4 and 3 Fr

**Length** 70,90,110,120 and 130 cm

