

# STONE BASKETS - SUPER CLENCH NITINOL RETRIEVER

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



## Key features:

- Multiwire design for tight clenching of fragments
- Partial closure forms a web of wires preventing retropulsion of stones
- Sturdy and kink resistant
- Unique design

## Specifications

### Size

**Diameter** 1.9,2.2,2.4 and 3 Fr

**Length** 70,90,110,120 and 130  
cm