


ULTRA HD EBUS Needle

An Endo-Bronchial Ultra Sound Needle system is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract. This device is supplied sterile and is disposable – intended for single use only. To be used by or on order of a physician only



Key features:

- Enhanced Echogenic needle tip for optimal ultrasonic visualization.
- Unique Needle design for higher and consistent tissue yield.
- Ergonomically contoured handle for better control and stability.
- Actuating lock for safer depth control of needle.
- Excellent needle flexibility for smooth tissue penetration.
- Ultra smooth sheath with adjustable length attuned for a full range of endoscopes.

Classification:

Minimum accessory channel 2.0 mm.

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com acomedicaldevices.com

Ref No.	Needle Description	Sheath Size
AEBUSN-19	19 gauge needle	4.1 FR
AEBUSN-21	21 gauge needle	4.1 FR
AEBUSN-22	22 gauge needle	4.1 FR
AEBUSN-25	25 gauge needle	4.1 FR

Specifications

Size

Sheath Diameter 4.1 Fr

Needle Size 19, 21, 22, 25

G