

# URETERAL STENTS OR DJ STENTS - CARBOTHANE

Used for internal drainage from the Uretropelvic junction to the bladder.  
Designed for mid-term indwelling.



## Key features:

- Highly radiopaque for better visualization
- Tapered at both the ends - for easy insertion
- Pigtail shape prevents stent migration
- Thin wall design for maximum drainage
- Smooth surface for easy insertion
- Available as a set with Guidewire, Clamp & Suture
- Available in White & Blue colour.

## Classification:

Models & Variants:

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email [info@acemedicaldevices.com](mailto:info@acemedicaldevices.com) [acemedicaldevices.com](http://acemedicaldevices.com)

- Open End / Closed End
- Both End Multiloop / One End Multiloop
- Stent with Guidewire
- Stent with Thread
- Stent with Clamp
- With / Without Hydrophilic Coating

## Specifications

### Size

**Diameter** 3 - 10  
Fr

**Length** 6 - 32 cm