




URETERAL STENTS OR DJ STENTS

Used for the temporary internal drainage from the ureteropelvic junction to the bladder.
Designed for short term use only.



Key features:

- Highly radiopaque for better visualization
- Tapered at both the ends - for easy insertion
- Pigtail shape prevents stent migration
- Thin wall design for maximum drainage
- Smooth surface for easy insertion
- Available as a set with Guidewire, Clamp & Suture
- Available in White & Blue colour.

Classification:

- Models & Variants:

Open End / Closed End

Both End Multiloop / One End Multiloop

Stent with Guidewire

Stent with Thread

Stent with Clamp

With / Without Hydrophilic Coating

Other properties:

Available in PolyUrethane, Carbothane Long Life and Silicon.

Specifications

Size

Diameter 3 - 10 Fr

Length 6 - 32 cm