



## URETERAL STENTS OR DJ STENTS

Used for the temporary internal drainage from the ureteropelvic junction to the bladder. Designed for short term use only.



### Key features:

- Highly radiopaque for better visualization
- Tapered at both the ends - for easy insertion
- Pigtail shape prevents stent migration
- Thin wall design for maximum drainage
- Smooth surface for easy insertion
- Available as a set with Guidewire, Clamp & Suture
- Available in White & Blue colour.

### Classification:

- Models & Variants:

Open End / Closed End

Both End Multiloop / One End Multiloop

Stent with Guidewire

Stent with Thread

Stent with Clamp

With / Without Hydrophilic Coating

## **Other properties:**

Available in PolyUrethane, Carbothane Long Life and Silicon.

## **Specifications**

### **Size**

**Diameter** 3 - 10 Fr

**Length** 6 - 32 cm