




Logo ace medical devices

URETERAL STENTS OR DJ STENTS

Used for the temporary internal drainage from the ureteropelvic junction to the bladder.
Designed for short term use only.



B1D35660 A063 444B 87E6C871937FCA9F

Key features:

- Highly radiopaque for better visualization
- Tapered at both the ends - for easy insertion
- Pigtail shape prevents stent migration
- Thin wall design for maximum drainage
- Smooth surface for easy insertion
- Available as a set with Guidewire, Clamp & Suture
- Available in White & Blue colour.

Classification:

- Models & Variants:

Open End / Closed End

Both End Multiloop / One End Multiloop

Stent with Guidewire

Stent with Thread

Stent with Clamp

With / Without Hydrophilic Coating

Other properties:

Available in PolyUrethane, Carbothane Long Life and Silicon.

Specifications

Size

Diameter 3 - 10 Fr

Length 6 - 32 cm