




URETHRAL STENTS (HYPOSPEDIAS STENT)

Used for stenting the urethra after hypospadias or Epispadias repair and to provide post operative drainage of the bladder. Designed for single use only.



Key features:

- Smooth surface being less traumatic
- Soft Polyurethane Material for patient comfort
- Highly radiopaque for better visualization
- Tapered end - for easy insertion.
- The loop design aids anchoring in bladder.
- Stiffening stylet for ease of insertion.

Specifications

Size

Diameter 6, 7, 8, 10, 12, 14, 16 Fr

Length 25 cm