3 STAGE BALLOON DILATOR OTW

Intended for endoscopic dilation of airway tree. Supplied sterile in a peel open pouch. Designed for single use only.



Key features:

- One balloon provides three different pre - defined diameters – 6-7-8mm, 8-9-10mm, 10-11-12mm, 12-13.5-15mm, 15-16.5-18mm
- Non compliant material offers optimum radial force and uniform inflated balloon diameter.
- Radiopaque bands to identify balloon location fluoroscopically
- Atraumatic tip prevents damage to the endoscope

- Accepts 0.035” guidewire

Specifications

Size

Sheath Diameter 6.8 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Length 1800 and 2400 mm

Balloon Length 55 mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com