
ANGIOGRAPHY GUIDE WIRE

Ace Angiography Guide Wires are hydrophilic guide wires and are used for percutaneous vessel entry to facilitate subsequent introduction(s) of an intravascular device. This device is Supplied sterile in a peel-open pouch. Intended for one-time use.



Key features:

- **Excellent trackability.**
- **Extremely smooth.**
- **Safe vessel access.**

Classification:

Straight and Angled Tip

Catalog No.

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

Size	Length	Type
ASFC 15018 H	0.018 Inches	150 cm
Straight		
ASFC 15021 H	0.021 Inches	150 cm
Straight		
ASFC 15025 H	0.025 Inches	150 cm
Straight		
ASFC 15032 H	0.032 Inches	150 cm
Straight		
ASFC 15035 H	0.035 Inches	150 cm
Straight		
ASFC 15038 H	0.038 Inches	150 cm
Straight		
AJFC 15018 H	0.018 Inches	150 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Angled Tip

AJFC 15021 H

0.021 Inches

150 cm

Angled Tip

AJFC 15025 H

0.025 Inches

150 cm

Angled Tip

AJFC 15032 H

0.032 Inches

150 cm

Angled Tip

AJFC 15035 H

0.035 Inches

150 cm

Angled Tip

AJFC 15038 H

0.038 Inches

150 cm

Angled Tip

Specifications

Size

Diameter 0.018 - 0.038 inch
G

Length 150 - 260 cm