
GUIDE WIRES - PTFE

AcePTFE Guidewires are used for percutaneous vessel entry to facilitate subsequent introduction(s) of an intravascular device. This device is Supplied sterile in a peel-open pouch. Intended for single use only..



Key features:

- Fixed core, Flexible floppy tip - for atraumatic advancement through urinory tract
- Stainless Steel with PTFE jacket - easy tracking of stents and uroligical instruments
- Excellent torqueability for precise tip control
- Low of coeffcient of friction for smooth advancement
- Better pushability for increased placement control
- Extremely smooth.
- Safe vessel access.

Classification:

Available in straight and J Tip variants

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

Specifications

Size

Diameter 0.018", 0.025", 0.028",0.032", 0.035" & 0.038"

Length 150 cm

Other

Balloon Max Fluid Injection

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com