ACE LAPAROSCOPIC SCISSORS DISPOSABLE

The Ace Indian Laparoscopic Monopolar Scissor is intended for use in Laparoscopic procedures involving cutting and coagulation.



Key features:

- High Carbon / Corrosion Resistance Stainless Steel Blades
- Dual- Edge Blades with a smooth cutting action
- Wide opening to take in larger chunks of tissue for cutting
- Does not loose its sharpness over extensive use
- Allows Dissection and transection of fine tissue
- High Quality Insulation TubingEnsuring insulation down to the blades (Dielectric Strength - 20 Kv/mm)
- Ergonomic Molded Handle that is comfortable for any size of hand.
- Molded with a Banana Monopolar Electrosurgery Connector
- Designed for a Comfortable 360 degree Rotation by your Fingertip.
- *India's FIRST and ONLY Manufacturer of Single Use Instruments.*

Specifications

Size

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

Diameter 5 mm

Working Length 330 ~ 430 mm

Jaw Length 18 mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com