
BILIARY METALLIC STENTS

Used for Stenosis or strictures in the Biliary tract. Supplied sterile in a peel open pouch with deployer.



Key features:

- Available as Fully Covered, Both Ends Bare, With removal string
- Available TTS and as Percutaneous transhepatic
- High corrosion resistance and shape memory
- Slow expansion below body temperature
- Unique deploying mechanism is reconstrainable over 98% of deployment to aid in accurate repositioning of the stent
- Optimum Radial Force
- Atraumatic stent endings to reduce tissue trauma
- Multiple radiopaque markers to facilitate accurate stent placement

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Specifications

Size

Diameter 6,8,10 mm

Length 40, 50, 60, 70, 80,100,120 mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com