BIOPSY GUN

Used for taking biopsies in CT, Ultrasound, surgery, urology, nephrology and other internal medicine applications.



Key features:

- Disposible Semi automatic Biposy needle
- Light weight and compact: spring loaded biopsy needle easier to manuever
- One hand operation:easily operated with one hand
- No risk of losing specimen: once the cannula is fired the Specimen is locked into the window of the stylet assuring an intact core sample for diagosis.
- Sharp beleved point allows quick easy penetration with less trauma to organs and surrounding tissues.

Classification:

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

14, 16 G x 10, 15 Cms - For Breast Biopsy
16, 18 G x 15 , 20 Cms - For Liver & Kidney
18, 20 G x 15, 20 Cms - For Lung
x 20, 25, 30 Cms - For Prostate

18 G

Specifications

Size

Diameter 14 - 18 G

Length 10 - 25 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com