CHIBA NEEDLE

Can be cosidered as a Universal Needle used for Percutaneous Nephrostomy procedures. Cytological biopsies needle aspirations thoracentesis exploratory punctures to inject contrast dyes.



Key features:

- Enables soft penetration and avoids the risk of vessel laceration and reaches deep points without any excessive deviation.
- The special echogenic tip increases the visibility under the ultrasound guidance and guarantees the needles point localization
- Bevel tip needle.

Specifications

Size

Diameter 18-22-23 G

Length 15 - 20 - 22 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com