# COLLINS KNIFE - BIPOLAR

Collin’s Knife is used to make straight incisions during TURP Surgery for treatment of BPH. Product supplied Sterile and single use only.



## Key features:

- Faster firing, longer lasting
- Excellent performance
- High quality insulation.
- Compatible with all standard resectoscopes

## Classification:

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

## Other properties:

Available with / without bipolar cable.

## Specifications

Size

**Diameter** 24
Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com