
CUTTING LOOP - BIPOLAR

Bipolar Cutting Loop is used to remove prostate tissue in the form of small chips during TURP surgery for treatment of BPH. Product supplied Sterile and single use only.



Key features:

- Faster firing, longer lasting.
- Excellent performance.
- High quality insulation.
- Compatible with all standard resectoscopes.

Classification:

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

Other properties:

Available with / without bipolar cable. Cutting Loops available in 0, 16, 20, 30, 45 and 90 Degrees.

Specifications

Size

Diameter 24
Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com