
ENUCLEATION BUTTON ELECTRODE - BIPOLAR

Ace Enucleation Button Electrode is used for bipolar enucleation, coagulation and vaporization for treatment of benign prostate hyperplasia (BPH). Compatible with Wolf and Storz 1 Stem Instrument. Product supplied Sterile and single use only.



Key features:

- Fast bipolar enucleation, coagulation and vaporization.
- Minimal penetration of thermal energy.
- Improved visibility and enhanced durability during procedures.

Classification:

Compatible with:

Storz 1 Stem Instrument

Wolf 1 Stem Instrument

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

Specifications

Size

Diameter 22, 24 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com