
FASCIAL DILATOR SET

Used for fascial dilation and intended for stretching or enlarging the fascial tissue covering a cavity, tract or opening prior to a diagnostic or interventional procedure. Designed for single use only



Key features:

- Smooth surface for ease of introduction
- Uniform taper to reduce trauma
- Can be supplied with or without hydrophilic coating

Specifications

Size

Diameter 6, 8, 10, 12, 14, 16, 18 and 20 Fr

Length 22 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com