FILIFORM URETHRAL DILATOR

Used for self dilation of Urethra.



Key features:

- Smooth surface for ease of introduction
- Uniform taper to reduce trauma.

Classification:

- With Hole
- Without Hole
- Curved

Specifications

Size

Diameter 6/6 - 22/6

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

Fr

Length 10 - 30 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com