
EUS - FNB NEEDLE

An Endoscopic Ultra Sound (EUS) - Fine Needle Biopsy (FNB) is used to sample targeted submucosal Gastrointestinal lesions through the accessory channel of an ultrasound endoscope. This device is supplied sterile and is intended for single use only. To be used by or on order of a physician only.



Key features:

- Enhanced Echogenic needle tip for optimal ultrasonic visualization.
- Unique needle design for higher and consistent tissue yield.
- Ergonomically contoured handle for better control and stability.
- Actuating lock for safer depth co ntrol of needle.
- Ultra smooth sheath with adjusta ble length attuned for a full range of endoscopes.

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com

Classification:

EUS Needle Ordering Information				
Needle Type	Needle Size Ø	19G	22G	25G
ACE SUPER CORE	FNB	AEUS-SFNB-19	AEUS-SFNB-22	AEUS-SFNB-25
ACE MASTER CORE		AEUS-MFNB-19	AEUS-MFNB-22	AEUS-MFNB-25
ACE ULTRA CORE		AEUS-UFNB-19	AEUS-UFNB-22	AEUS-UFNB-25
ACE DUAL CORE	FNB	AEUS-DFNB-19	AEUS-DFNB-22	AEUS-DFNB-25
ACE SUPER DUAL CORE		AEUS-SDFNB-19	AEUS-SDFNB-22	AEUS-SDFNB-25
ACE MASTER DUAL CORE		AEUS-MDFNB-19	AEUS-MDFNB-22	AEUS-MDFNB-25
ACE ULTRA DUAL CORE		AEUS-UDFNB-19	AEUS-UDFNB-22	AEUS-UDFNB-25

Specifications

Size

Sheath Diameter 5.2 Fr

Needle Size 19, 22 and 25G

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com