
HYDROPHILIC GASTRO GUIDEWIRE

Used to facilitate the placement of instruments during diagnostic or treatment procedures. Supplied sterile in a peel open Tyvek pouch. For single use only



Key features:

- Unmatched torqueability
- Excellent trackability.
- Extremely smooth
- Safe vessel access
- Nitinol non kink core

Specifications

Size

Diameter 0.018, 0.021, 0.025, 0.028, 0.032, 0.035 Inch

Length 260 and 450 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com