INITIAL PUNCTURE NEEDLE 2 PART

Used for Percutaneous Nephrostomy, cyst puncture, deep tissue access and all other Percutaneous application. Designed for single use only.



Key features:

- Diamond & Bevel tip

Specifications

Size

Diameter 18 G

Length 15 - 23 Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com, acemedicaldevices.com