# MECHANICAL LITHOTRIPTOR

Used for regular or emergency Lithotripsy to crush C.B.D. stones endoscopically. Reusable device supplied non sterile.



## Key features:

- Available as TTS and OTW variants
- Spirals ensure smooth and easy lithotripsy
- Basket shape hexagonal
- Spiral OTW : 4 mm, length 100 cms
- Spiral TTS : 3.5 mm, length 195 cms
- Reusable design that can be steam sterilized
- Supplied in a carrying case

## Specifications

Size

**Basket Diameter** 25, 30 and 35 mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com

**Length** 400 cm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com, acomedicaldevices.com